

# Application Note



Akademie věd České republiky  
Ústav teorie informace a automatizace AV ČR, v.v.i.

## Yeast Colony Counter User Manual

Jan Schier, Bohumil Kovář  
[schier@utia.cas.cz](mailto:schier@utia.cas.cz), +420-2-6605 2470

*Centre for Applied Cybernetics  
1M0567*

*Department of Signal Processing  
<http://zs.utia.cz>*

### Contents

<b>1 Introduction</b>	<b>1</b>
1.1 Image characteristics . . . . .	1
<b>2 Installation of the tool</b>	<b>1</b>
<b>3 Launching the tool</b>	<b>1</b>
<b>4 Home screen</b>	<b>2</b>
4.1 First start of the tool . . . . .	2
4.2 Subsequent run of the tool . . . . .	3
4.3 Selection of the result file . . . . .	3
<b>5 Processing the dishes in the supervised mode</b>	<b>4</b>
<b>6 Result file – filename and format</b>	<b>5</b>
<b>7 Processing flow</b>	<b>5</b>
<b>8 Image checking</b>	<b>6</b>
<b>9 A note on the tool performance</b>	<b>6</b>
<b>10 Acknowledgements</b>	<b>6</b>

## Revize

Revize	Datum	Autor	Popis změn v dokumentu
0	1.10.2010	J.S.	Vytvoření dokumentu

# 1 Introduction

This document describes the use of the yeast colony counter developed in UTIA AV CR, Prague. This colony counter provides batch processing of a series of images of Petri dishes containing yeast colonies. It is aimed for a research laboratory using standard imaging equipment (darkroom, general purpose camera on a stand) and hand manipulation of the dishes rather than specialized single-purpose camera system. It is assumed that dark background is used in the images.

The tool evaluates two figures for the images: the relative coverage of the dish, and the number of colonies in the dish. It provides two modes of operation: automatic or semi-automatic. In the automatic mode, it processes the batch of images without any user intervention. In the semi-automatic mode, it first performs automated thresholding, dish localization, and a counting attempt, then it provides the user with an editor, where he can interactively correct the result of counting.

Two image formats are accepted at this moment: jpg and tiff. Other formats may be added on request.

## 1.1 Image characteristics

The images under consideration are characterized by the following properties:

- + typical resolution of an image is  $1024 \times 768$  pixels,
- + dark background is used to increase image contrast,
- + illumination by two linear lights along the short edges,
- + colonies are roughly round-shaped,
- variations of dish position,
- varying diameter of the dish in image,
- varying intensity of the dish background.
- colonies often form clusters.

## 2 Installation of the tool

The tool comes in the form of a self-extracting archive, which contains the tool self (ColonyCounter.exe), and the Matlab Compiler Runtime Installer (MCRInstaller.exe), which is necessary to run the tool. When running the ColonyCounter\_pkg archive, the files are extracted into the directory of archive and the MCRInstaller is automatically launched. It asks to choose the setup language (choose the default – “English”) and for permission to install the VCREDIST\_X86 (Visual C++ Redistributable) (choose “Install”).

The installation continues automatically until the MCR InstallShield Wizard asks for several installation options: fill in the User information or just click on “Next”; select installation directory and then launch the installation. The InstallShield may ask about installation of the .NET Framework - just click “OK” to continue.

When the installation process is finished, the *ColonyCounter.exe* program is unpacked to the directory. This is the tool self — other files can be safely removed.

## 3 Launching the tool

To launch the tool, simply run the ColonyCounter executable, e.g. with *Start*→*Run*→*ColonyCounter*. The start of the program takes some 15 s on our system.

## 4 Home screen

After launching the tool, the home screen appears, as shown in Figure 1. It is used to select the following:




Figure 1: Home screen of the tool

- *Image directory*: directory with the images to process (directory tree is possible as well)
- *Result directory*: location and file name for the file containing the results of processing
- *Mode*: either *Supervised* — with graphical editor for correction of the colony detection results, or *Batch* — fully automatic.
- *Start processing*: is used to start the processing, after making the necessary selections.

### 4.1 First start of the tool

On the first start of the tool, the *Result directory* and *Start processing* items are disabled. The user is forced to make the selection of the image directory, after which the *Result directory* item will be enabled. Finally, after having selected the result file, the *Start processing* item will also get enabled.

**Note:** the values of the image directory and result file selections are stored in the *ColCounterConfig.txt* file between the runs of the tool.

## 4.2 Subsequent run of the tool

On the subsequent runs of the tool, the values of the image directory and result file from the last run are re-loaded. The *Start processing* item is now enabled on the tool launch. If the image directory is changes, warning about the change is issued, reminding the user to check the location of the result file (see Figure 2), and the *Start processing* button gets disabled.




Figure 2: Warning message after directory change

## 4.3 Selection of the result file

The default mode, which is used if the user does not make any selection on a subsequent run of the tool (neither image directory nor result file), is to *continue* the processing of the given image directory and to *append* to the result file. The tool goes through the result file, checks the results against the files in the directory tree, delete the data for the images that are not present any more, and continues with processing at the point where is has been interrupted on the last run.

If the user decides to open the *Result file* selection window, and *selects the same file*, he gets a series of *warning dialogues*:

- First warning about overwriting the file comes from MSWindows and could not be eliminated. Nevertheless, the user will have chance to *choose to append or overwrite* the results file in the following dialogue windows, so it is safe to choose to overwrite at this point.
- The warning that is issued next, is shown in Figure 3. Here, user is offered the possibility either to




Figure 3: Warning dialogue: existing result file

*Append* (default) or to *Overwrite* the result file.

- If he selects to overwrite the file, he is warned once more, as shown in Figure 4, with the default selection set to *No* (that is, do not overwrite the results). Selecting *Yes* in this dialogue will result in overwriting the results.




Figure 4: Overwrite conformation dialogue

## 5 Processing the dishes in the supervised mode

In the supervised mode, the results of processing (the detected colonies) are presented in an editor window, as shown in Figure 5. This window has several function and information areas:




Figure 5: Dish editor window

- *Window title*: path and filename of the image
- *Top area*: results: relative area and number of colonies
- *Central area*: image of the dish with colonies, with marks on detected colonies
- *Bottom part*: Edit buttons *Delete marks*, *Add marks*, *Zoom* and *Next dish*.

The meaning of *Delete* and *Add* is obvious. *Zoom* button switches to zoom mode, which can be used to zoom-in the area-of-interest in the dish. Both click-hold-drag and scroll-wheel can be used. *Next dish* is used to finish processing of the given dish and to proceed to the next one.

- *Right column*: *Stop* button can be used to interrupt processing. *Save image* can be used to save image of the dish with detection signs for archival purposes.

**Note:** the filename of the saved images is, by default, given an appendix `_saved`. The user can choose any filename and location; however, the images with this suffix, if present in the image directory, are ignored by the tool when it seeks for the images to process (obviously, processing screenshots would give little sense).

## 6 Result file – filename and format

The file with results is standard CSV (Comma-Separated Values) file, which can be loaded into any spreadsheet program. The values are separated by semicolon.

The default name of this file, as suggested by the tool, is the name of the root directory of the image tree with the `_results` suffix, e.g. `images_results.csv`. However, any filename can be selected by user.

A line in the file has the following structure:

1. *path to the image* from the root of the image directory tree (e.g. `images/dense/1.jpg`). Note: this is the same value as displayed in the title-bar of the image editor (Sec. 5).
2. *relative coverage* of the dish in percents
3. *number of colonies* in the dish

The result file is updated immediately after processing each dish, so it is save to interrupt the processing.

## 7 Processing flow

Processing flow of the tool is outlined in Figure 6.




Figure 6: Processing flow of the colony counter

First step of the processing includes thresholding of the image and localisation of the dish — let us recall that the dishes are hand manipulated (Sec. 1), hence, the position of the dish in each image is slightly different. In this step, the image is also checked for correctness.

The thresholding is used to separate the dish from the background and to perform image segmentation, that is, to extract the (clusters of) blobs.

Next step is resolution of the clustered blobs: the Fast Radial Transform [Loy and Zelinsky, 2003] is used to construct symmetry map of the objects in the image. Then, the nonmaxsuppts function [Kovesi, 2005] is used to find local maxima of this map, which represent the centers of the colonies.

Finally, the position of the centers is checked to filter out cases which would be obviously wrong (e.g. colony center falling in the background).

## 8 Image checking

In the first step, image is also checked for correctness, to ensure sufficient robustness of the processing.

This check involves mainly the position and size of the dish: it should be roughly centered and the inner part (that is, the part with the grow medium) must be completely present in the image, so that the coverage is evaluated correctly.

The image of the dish must be also sufficiently large, to provide reasonable accuracy of detection. Currently, the setting is that the dish must occupy at least 50% of the image height.

The background distortions, such as that shown in Figure 7, are eliminated in the process of dish location and background masking, as long as the distortion does not overlap with the dish in vertical projection.



Figure 7: Example of background error

## 9 A note on the tool performance

The performance of the tool has been tested using 245 images of Petri dishes, containing mostly less than 60 colonies per dish. The performance is summarized in Table 1. More details are given in a conference paper that is currently under review process. A reference will be supplied as soon as possible.

## 10 Acknowledgements

The images included in the package have been prepared in the Department of Genetics and Microbiology, Faculty of Natural Sciences, Charles University. The software has been developed in cooperation with this laboratory.



# colonies	samples	missed [%]	false [%]
0 – 17	21	1.47	0
18 – 20	19	2.99	0
21 – 23	25	3.31	0
24 – 26	29	4.00	0
27 – 29	27	3.28	0
30 – 40	35	4.17	0
41 – 49	30	4.00	0
50 – 60	29	3.87	2.00
> 60	30	5.36	3.22

Table 1: Dependence of the counting errors on the number of colonies in the dish. The system average counting error is under 4%

## References

- [Kovesi, 2005] Kovesi, P. D. (2005). MATLAB and Octave functions for computer vision and image processing. School of Computer Science & Software Engineering, The University of Western Australia. Available from: <http://www.csse.uwa.edu.au/~pk/research/matlabfns/>.
- [Loy and Zelinsky, 2003] Loy, G. and Zelinsky, A. (2003). Fast radial symmetry for detecting points of interest. *IEEE Transactions on Pattern Analysis and Machine Intelligence*, 8(25):959–973.